## **Wood Carbonization Protocol**

This protocol is for the carbonization of modern wood specimen, to accession into the EAL and JMM wood charcoal comparative collections. It will take approximately 10-40 minutes to carbonize specimen in the muffle furnace, but allocate 1-2 hrs to wrap, carbonize, and unwrap 10 specimens.

#### SAMPLE PREPARATION

1. Be sure the samples are accessioned into the FileMaker database and assigned an EAL accession number.

#### PREPARING SPECIMEN ~ 20 - 60 mins

- 1. Cut or saw wood into pieces so that the smallest dimension is no more than 1.5-2.0 inches thick.
- 2. Using a pencil, write the accession number in multiple locations on the specimen.
- 3. Wrap in aluminum foil as tightly as possible, leaving a small flap/tag of foil on to which you will write the accession number.
- 4. Using an Ultra Fine Sharpie, engrave the accession number onto the flap and onto the foil. The ink will burn away during firing but the engraved number will remain.

#### CARBONIZING SPECIMEN ~ 15 - 45 mins

- 1. Turn on the fume hood.
- 2. Place specimen inside the muffle furnace, with larger specimen on the bottom and smaller ones closer to the front.
- 3. Draw a diagram of specimen locations, to ensure proper identification following firing.
- 4. Be sure the ceramic stopper is in the top of the muffle furnace to reduce oxygenated atmosphere.
- 5. Set muffle furnace to 400°C (725°F).
- 6. Fire specimens between 10-40 mins.
- \*Be sure not to raise the sash higher than approved by Health and Safety when working inside fume hood and close sash fully when carbonizing.

### REMOVING SPECIMEN ~ 15 mins

- 1. Remove carbonized specimen when cool, and examine to see if complete carbonization was achieved.
- \*If specimens are still hot, an introduction of oxygen could ignite a fire during unwrapping. Ensure specimens are cool before unwrapping.
- 2. Turn off muffle furnace and fume hood.
- 3. Unwrap tin foil and place specimens in labeled individual cardboard sample boxes with their corresponding accession number.

# **Revision History:**

Created by Kathleen M. Forste, Boston University: 1/30/17 Standardized by Kali Wade, Boston University: 1/17/18