Boston University Sargent College of Health & Rehabilitation Sciences Speech, Language & Hearing Sciences

BACKGROUND

- Speech understanding in multitalker mixtures requires both segregation of the target from distractors and identification based on partial information (or "speech glimpses").
- Identification based on speech glimpses is particularly susceptible to subtle reductions in bandwidth (Best et al 2019) and audibility (Best et al 2017).
- An implication is that restoring audibility across the spectrum for individuals with hearing loss is especially critical for identification of speech in multitalker mixtures.

AIMS & RATIONALE

- Here we asked whether current hearing-aid amplification strategies adequately restore the audibility of speech glimpses.
- In a multitalker mixture, gain and compression may often be driven by the maskers and it is difficult to know how audible the target actually is.
- There may also be opposing effects that counteract improvements in target audibility (Overby et al 2023):
 - amplification of maskers may make them more distracting and harder to ignore
 - non-linear processing may distort natural cues for segregation of competing talkers (e.g., cross-modulation, spatial cue distortion).
- Are these opposing effects why hearing aids are often not helpful in real-world mulitalker situations?

REFERENCES

Best et al (2019). The importance of a broad bandwidth for understanding "glimpsed" speech. J Acoust Soc Am 146:3215-3221.

Best et al (2017). Use of a glimpsing model to understand the performance of listeners with and without hearing loss in spatialized speech mixtures. J Acoust Soc Am 141:81-91.

Overby et al (2023). A data-driven distance metric for evaluating the effects of dynamic range compression in adverse conditions. Trends Hear 27:23312165231192302.

Brungart et al (2006). Isolating the energetic component of speech-on-speech masking with ideal time-frequency segregation. J Acoust Soc Am 120:4007-4018.

Kayser et al (2022). Open community platform for hearing aid algorithm research: open Master Hearing Aid (openMHA). SoftwareX 17:100953.

Research reported in this publication was supported by the National Institute On Deafness And Other Communication Disorders of the National Institutes of Health under Award Number R01DC015760. The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

How well do hearing aids amplify "speech glimpses" in multitalker mixtures?

Virginia Best (ginbest@bu.edu)





used to isolate speech glimpses from unaided and

aided speech mixtures.

- dB SPL) and two masker sentences (male, variable level) presented via headphones.
- Individualized non-linear amplification (NAL-NL2) applied using a hearing-aid simulator (Oldenburg openMHA; Kayser et al 2022).
- Participants were 10 younger adults (age 19-40) with bilateral sensorineural hearing loss. 7 were regular hearing-aid users.

- On average, amplification did improve the audibility and intelligibility of the target speech glimpses. These improvements were reduced in the mixture, suggesting that there were counteracting effects.
- On an individual level, this pattern was true for about half of the participants. For others, amplification did not clearly improve the intelligibility of target speech glimpses.
- This approach provides a useful way to unpack positive and negative effects of hearing-aid processing in multitalker mixtures.



APPROACH

Ideal time-frequency separation (Brungart et al 2006)

Stimuli were mixtures of a target sentence (female, 65)



RESULTS

-Mix Aided Glimpsed Unaided Glimpsed Aided Quiet Target-to-masker ratio (dB)



